

# Steris Amsco 3080 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Steris Amsco 3080 Premium Series Replacement Pad Set is a 3-piece cushion set engineered for durability and consistent support during long surgical procedures. It fits Steris/Amsco surgical table models including the 3080RC, 3080RL, 3080SP, and 3085SP.

This replacement pad set balances support and comfort using high-quality materials. The Visco memory foam layer conforms to the patient's body, reducing pressure points. The Fusion III vinyl cover stretches in four directions, improving durability and flexibility while its conductive surface helps prevent static buildup in the O.R. environment.

The combination of medium-firm high-density foam with memory foam and premium vinyl covering makes this pad set ideal for hospitals that demand long-term performance and reliable patient positioning during demanding operations.

## FEATURES

- For Longer Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)