

Steris Amsco 3085 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Steris Amsco 3085 Standard Series Replacement Pad Set is a 3-piece cushion system designed for short procedures and low-risk patients. It fits Steris/Amsco models 3080RC, 3080RL, 3080SP, and 3085SP, providing an affordable solution for surgical teams requiring basic patient support.

The pad set features a durable high-density foam core that maintains patient stability during short-duration cases. The 2-inch thick cushions are constructed with medium-firm foam for reliable support during routine procedures.

The electrically conductive vinyl cover resists static, fluids, and disinfectants, making it suitable for everyday surgical environments. The lightweight design and easy installation help reduce setup time between procedures.

FEATURES

- For Low-Risk Patients
- For Short Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Hook & Loop Kit Included

SPECIFICATIONS

- Construction: Medium Firm High-Density Exterior
- Cover: Electrically Conductive Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)