

Steris Amsco 3085 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Steris Amsco 3085 Premium Series Replacement Pad Set is a 3-piece pad set designed for longer surgical procedures. It fits Steris/Amsco models 3080RC, 3080RL, 3080SP, and 3085SP, and can be used to replace worn pads and restore surgical table performance.

Each pad combines medium-firm high-density foam with 1 inch of Visco memory foam, totaling 2 inches in thickness. The Visco memory foam layer adapts to body contours for enhanced comfort and pressure relief. The pads are covered with Conductive 4-way stretch Fusion III Soft Vinyl that resists fluid and stains while withstanding hospital-grade cleaners.

The conductive surface reduces static in the operating room environment. The pad set enhances patient support while reducing pressure-related risks during extended procedures.

FEATURES

- For longer Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)