

Steris 4085 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Steris 4085 Surgical Table Standard Series Replacement Pad Set is a 3-piece cushion set designed for low-risk patients and short-duration procedures. It fits both Steris 4085 and CMax surgical tables, providing a versatile and cost-effective solution for general operating room use.

The pad set features 2-inch thick medium-firm high-density foam construction covered with electrically conductive vinyl. This design delivers consistent support for routine surgical procedures while the conductive vinyl cover helps reduce static and resists fluids and disinfectants commonly used in hospital environments.

The lightweight design ensures easy installation and daily usability in fast-paced operating room settings. The pads are constructed to provide reliable comfort and support while maintaining durability in demanding clinical environments.

FEATURES

- For Low-Risk Patients
- For Short-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free
- Hook and Loop Kit included

SPECIFICATIONS

- Construction: Medium Firm High-Density
- Exterior Cover: Electrically Conductive Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)