

# Steris 4085 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Steris 4085 Premium Series Replacement Pad Set is a 3-piece cushion system designed for extended surgical procedures. It fits both Steris 4085 and Steris CMax surgical tables, providing versatile compatibility across operating room environments.

The pad set features 2" thick construction combining medium-firm high-density foam with 1" of visco memory foam to relieve pressure and enhance patient positioning during long procedures. The conductive 4-way stretch Fusion III Soft Vinyl cover resists fluid and stains while withstanding hospital-grade disinfectants.

The conductive surface reduces static buildup to support a safer operating environment for both staff and patients. This premium pad set prioritizes patient safety and comfort while maintaining durability in demanding surgical settings.

## FEATURES

- For Longer-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)