

# Steris CMax Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Steris CMax Standard Series Replacement Pad Set is a 3-piece cushion set designed for short-duration procedures and low-risk patients. This pad set fits both the Steris CMax and Steris 4085 surgical tables, providing flexibility for multi-table facilities. The high-density foam core maintains patient positioning throughout the procedure while the conductive vinyl cover reduces static and resists fluids and harsh disinfectants.

Constructed with medium-firm high-density foam and covered with electrically conductive vinyl, this pad set provides consistent support and stability in general surgical settings. The lightweight design and easy installation make it practical for everyday operating room use. The conductive vinyl cover is fluid and stain resistant, helping to protect against common surgical environment challenges.

## FEATURES

- For Low-Risk Patients For Short-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Hook & Loop Kit included

## SPECIFICATIONS

- Construction: Medium Firm High-Density
- Exterior Cover: Electrically Conductive Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)