

# Steris CMax Premium Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Steris CMax Premium Series Replacement Pad Set is a 3-piece cushion system designed to improve patient positioning and pressure relief during extended surgical procedures. It fits both the Steris CMax and Steris 4085 surgical tables. The set features 2" thick construction with medium-firm high-density foam and 1" visco memory foam to reduce pressure points and enhance patient stability.

The pads are covered with Conductive 4-way stretch Fusion III Soft Vinyl that resists fluid, tearing, and disinfectants. The conductive surface helps reduce static electricity, providing an added layer of safety in surgical environments. The cover stretches in four directions for optimal fit and comfort during positioning.

## FEATURES

- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)