

Steris 5085 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Steris 5085 Premium Series Replacement Pad Set is a 3-piece cushion system designed for extended surgical procedures in high-demand operating room environments. The set features 2 inches of medium-firm high-density foam with 1 inch of Visco memory foam, topped with a conductive 4-way stretch Fusion III Soft Vinyl cover. This construction enhances pressure management, improves patient positioning, and helps reduce the risk of pressure injuries during longer operations.

The Visco memory foam top layer conforms to the patient's body, increasing comfort and stability during procedures. The Fusion III Soft Vinyl cover provides four-way stretch for flexibility and durability while resisting fluid penetration, tearing, and hospital disinfectants. The conductive surface reduces static buildup, adding an important safety layer in sensitive surgical settings.

This pad set is built for durability and reliable performance in demanding surgical environments. It delivers consistent support and comfort across extended procedures while maintaining the structural integrity required for professional medical use.

FEATURES

- For Longer-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)