

# Steris 5085 SRT Premium Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Steris 5085 SRT Premium Series Replacement Pad Set is a 3-piece cushion system designed to provide reliable support, pressure relief, and positioning during surgical procedures. The set combines 1" of visco memory foam with a medium-firm high-density foam core, totaling 2" in thickness, to deliver lasting comfort and support throughout extended operations.

The pads are covered with Conductive 4-way stretch Fusion III Soft Vinyl, which offers superior stretch, fluid resistance, and durability. This cover withstands daily cleaning with hospital-grade disinfectants and reduces static buildup in the operating room environment.

The pad set is suitable for hospitals and surgical centers seeking to maintain high performance and safety standards. The combination of materials and construction provides effective pressure relief while meeting the demands of long and challenging surgical procedures.

## FEATURES

- For Longer-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)