

Steris SurgiGraphic 6000 Standard Series Replacement Pad

FACTORY NEW

The Steris SurgiGraphic 6000 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. This cost-effective option provides essential support and surface protection during routine imaging and surgical cases. Hospitals and outpatient centers use this pad set when extended pressure relief is not required.

The pad set features a 2-inch thick core constructed with medium-firm high-density foam that maintains shape under light to moderate use. The electrically conductive vinyl cover helps reduce static buildup in imaging environments and resists common cleaning agents.

The lightweight design allows for easy installation and handling. The high-density foam core delivers reliable stability and basic patient comfort for routine surgical and imaging applications.



Stock Photo - Options May Vary

FEATURES

- For Low-Risk Patients For Short Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Hook & Loop Kit Included

SPECIFICATIONS

- Construction: Medium Firm High-Density
- Exterior Cover: Electrically Conductive Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)