

# Steris SurgiGraphic 6000 Premium Series Replacement Pad

## FACTORY NEW

The Steris SurgiGraphic 6000 Premium Series Replacement Pad Set is engineered for patient comfort during longer-duration surgical procedures. This high-quality cushion set offers reliable pressure relief and support throughout extended operations.

The pad features a 2" thick construction with medium-firm high-density foam and 1" Visco memory foam layer that molds to the patient's shape, improving comfort and reducing shear forces. The Visco memory foam layer molds to the patient's shape, improving comfort and reducing shear forces.

The Fusion III Soft Vinyl cover stretches in four directions, offering durability and flexibility while resisting fluids, tears, and harsh disinfectants. The conductive surface minimizes static buildup for added safety in the operating room, making it ideal for busy surgical settings.



Stock Photo - Options May Vary

## FEATURES

- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)