

MCM Bariatric Drop-latch Extra Wide Armboard | MCM 412



Stock Photo - Options May Vary

FACTORY NEW

The MCM 412 Drop-Latch Style Bariatric Armboard is designed for strength and reliability in surgical settings. Setup is fast—simply drop the board onto the rail where a weighted trigger lock instantly secures it in place for safe and stable support. Removal is equally straightforward: hold the trigger lock up and lift the outer end of the board to release.

The finger latch trigger, located at the outer end, allows for a full 180-degree range of rotation, giving surgical teams better control and positioning during procedures. Made from durable 3/8-inch thick black phenolic with a rugged steel gear casting, this armboard offers excellent strength without unnecessary weight.

Measuring 10 inches wide and 26 inches long, the board fits all standard American-made OR tables and accommodates larger patients. The armboard pad is sold separately.

FEATURES

- Armboard drops onto rail and locks by means of a weighted trigger lock.
- Armboard releases off rail by simply holding trigger lock up and lifting on the outer end of board.
- Finger latch trigger on outer end of board allowing for 180 degree rotation
- Durable steel casting for longevity
- Armboard is made of 3/8" thick black phenolic and weighs a slim 4.5 pounds
- Pad sold separately (MCM 413)
- 10" Wide x 26"L

SPECIFICATIONS

- Length: 26"
- Width: 10"
- Made of 3/8" thick black phenolic
- Weight: 6.5 pounds
- Dimensions: 10" W x 26" L

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)