

# MCM Standard Drop-Latch Armboard – 24"L | MCM 410



Stock Photo – Options May Vary

## FACTORY NEW

The MCM 410 is a standard drop-latch surgical table armboard manufactured by Mid Central Medical for use in operating rooms and imaging environments. It features a convenient drop-on rail design with a weighted trigger mechanism that allows clinicians to quickly secure and release the armboard as needed during procedures.

The armboard incorporates a finger latch on its outer edge that enables smooth 180° lateral rotation for accurate arm positioning. Its reinforced gear casting eliminates the risk of broken teeth, ensuring long-lasting structural integrity and reliable performance in demanding surgical settings.

Constructed from 3/8" thick black phenolic material, the MCM 410 is compatible with all standard American-made operating room tables. At just 4.5 pounds, it provides a lightweight yet durable solution for surgical positioning needs.

## FEATURES

- Armboard drops onto rail and locks by means of a weighted trigger lock
- Armboard releases off the rail by simply holding the trigger lock up and lifting the outer end of the board
- Finger latch trigger on the outer end of the board allowing for 180° lateral adjustment
- Designed with extra strong gear casting to eliminate broken teeth on castings
- Armboard is made of 3/8" thick black phenolic
- Weighs a slim 4.5 pounds
- Fits all standard American made operating room tables
- Tried and proven armboard

## SPECIFICATIONS

- Armboard Length: 24"
- Armboard Width: 6"
- Armboard Weight: 4.5 pounds

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)