

# MCM Standard Drop-Latch Armboard - 26"L | MCM 411



*Stock Photo - Options May Vary*

## FACTORY NEW

The MCM 411 is a 26-inch standard surgical armboard designed for use with North American surgical table side rails. It features a drop-latch design with a weighted trigger lock for quick attachment and removal without tools. The board supports 180-degree lateral rotation via a finger latch trigger on the outer end, enabling fast and controlled positioning during procedures.

Constructed from durable 3/8-inch thick black phenolic material, the armboard weighs only 4.5 pounds while maintaining strength and stability. The enhanced locking mechanism ensures firm positioning during use, and hook and loop material is pre-attached to secure the pad in place.

The armboard is built for reliability in daily surgical use with a strong gear casting designed to prevent damage. It offers universal compatibility with standard surgical tables and comes with a two-year warranty.

## FEATURES

- Armboard drops onto rail and locks by means of a weighted trigger lock.
- Armboard releases off rail by simply holding trigger lock up and lifting on the outer end of board.
- Finger latch trigger on outer end of board allowing for 180 degree rotation
- Durable steel casting for longevity
- Armboard is made of 3/8" thick black phenolic and weighs a slim 4.5 pounds
- Pad sold separately
- Available in 24" or 26" long sizes.

## SPECIFICATIONS

- Length: 26"
- Width: 6"
- Made of 3/8" thick black phenolic
- Weight: 4.5 pounds
- Dimensions: 6" W x 26" L

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)