

MCM Standard Drop-Latch Armboard with 2" Pad - 26"L | MCM 411-W2P



Stock Photo - Options May Vary

FACTORY NEW

The MCM 411-W2P is a 26-inch armboard designed for surgical tables with standard North American side rails. It features a drop-latch design with a weighted trigger lock for quick attachment and removal without tools. The board rotates 180 degrees laterally via a finger latch trigger on the outer end, enabling fast and controlled positioning.

Constructed from 3/8-inch thick black phenolic material, the armboard weighs only 4.5 pounds while maintaining durability and strength. Its enhanced locking mechanism ensures firm and consistent positioning during surgical procedures. Hook and loop material is pre-attached to secure the included 2-inch pad.

The armboard is built for everyday operating room use with a universal fit design and durable steel casting components. Its lightweight construction and simple setup make it an efficient addition to any surgical table system.

FEATURES

- Armboard drops onto rail and locks by means of a weighted trigger lock.
- Armboard releases off rail by simply holding trigger lock up and lifting on the outer end of board.
- Finger latch trigger on outer end of board allowing for 180 degree rotation
- Durable steel casting for longevity
- Armboard is made of 3/8" thick black phenolic and weighs a slim 4.5 pounds
- 2" Pad comes included.

SPECIFICATIONS

- Length: 26"
- Width: 6"
- Made of 3/8" thick black phenolic
- Weight: 4.5 pounds
- Dimensions: 6" W x 26" L

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)