

# MCM Standard Drop-Latch Armboard with 2" Pad - 24"L | MCM 410-W2P



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 410-W2P is a 24" surgical armboard designed for standard North American operating table rails. It features a drop-latch weighted trigger mechanism for quick, tool-free installation and removal. The board supports 180-degree lateral rotation via a finger latch trigger on the outer end, enabling precise positioning in any direction.

Constructed from durable 3/8-inch thick black phenolic material, the armboard weighs only 4.5 pounds while maintaining strength and stability. The enhanced locking mechanism ensures firm, consistent positioning during surgical procedures. A strong gear casting helps prevent tooth breakage and extends the board's service life.

The armboard measures 26 inches long and 6 inches wide and includes a 2-inch pad with pre-attached hook and loop material. It is compatible with standard surgical tables and backed by a two-year warranty.

## FEATURES

- Armboard drops onto rail and locks by means of a weighted trigger lock.
- Armboard releases off rail by simply holding trigger lock up and lifting on the outer end of board.
- Finger latch trigger on outer end of board allowing for 180 degree rotation
- Durable steel casting for longevity
- Armboard is made of 3/8" thick black phenolic and weighs a slim 4.5 pounds
- 2" Pad is Included.

## SPECIFICATIONS

- Length: 24"
- Width: 6"
- Made of 3/8" thick black phenolic
- Weight: 4.5 pounds
- Dimensions: 6" W x 24" L

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)