

# MCM Extra Wide Underpad Armboard | MCM 442



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The MCM 442 is an extra wide underpad armboard designed to provide enhanced support and positioning for surgical and imaging procedures. Its 21" x 13.25" surface accommodates various patient sizes while offering improved comfort and stability during procedures.

Constructed from durable, radiolucent material, the MCM 442 is ideal for C-arm imaging and fluoroscopic procedures. The radiolucent design ensures clear, unobstructed imaging, making it suitable for operating rooms and interventional suites where precision and visibility are critical.

The armboard slides under the table pad and remains in place using the patient's weight, requiring no clamps or fasteners. Its low, simple design provides steady support without interfering with other equipment, making it a reliable solution for both basic positioning and detailed imaging applications.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)