

MCM UnderPad Mount Armboard - 24" L | MCM 440



Stock Photo - Options May Vary

FACTORY NEW

The MCM 440 Underpad Mount Armboard is designed for tables without standard side rails or where imaging access is critical. It slides beneath the table pad and is held in place by the patient's weight, requiring no clamps or attachments. This design enables infinite positioning options, giving surgical teams total freedom to adjust for comfort and access.

Constructed from durable 3/8-inch thick phenolic, the armboard maintains rigidity and strength without needing extra support. It weighs just 7 pounds, making it easy for nurses and staff to position it wherever needed.

The armboard is fully radiolucent with no metal parts to interfere with imaging equipment such as C-Arms. A compatible pad is sold separately and should match the thickness of your existing table surface.

FEATURES

- Armboard is made of 3/8" thick black phenolic and weighs a mere 7 pounds.
- Fully Radiolucent

SPECIFICATIONS

- Length: 24"
- Width: 6"
- Made of 3/8" thick black phenolic
- Weight: 7 pounds
- Dimensions: 6" W x 24" L

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)