

# MCM Double Picket Fence Leg Holder | MCM 610

## FACTORY NEW

The MCM 610 Double Picket Fence Leg Holder is designed for secure and efficient leg positioning during surgical prep. Commonly used in orthopedic and vascular procedures, it features a durable urethane-covered top for added patient comfort and easy cleaning. The unit attaches with two MCM 162 flat blade sockets (sold separately), offering a secure and adjustable setup.



*Stock Photo - Options May Vary*

This versatile leg holder fits any fluoroscopy or imaging table, making it ideal for various clinical environments. It supports arms, legs, and feet in multiple positions, helping improve procedural accuracy and efficiency. The double picket fence design enhances stability while allowing for quick adjustments during patient preparation.

Its compact 7.5-inch wide frame conveniently folds for easy storage, reducing clutter in the operating room. Built to last and simple to use, the MCM 610 is a reliable solution for modern surgical suites needing adaptable, space-saving patient support equipment.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)