

MCM Stainless Steel Clark Socket | MCM 160



Stock Photo - Options May Vary

FACTORY NEW

The Mid Central Medical MCM 160 is a stainless steel Clark Socket designed for use in surgical environments. It accepts round posts and offers versatile attachment options for surgical accessories. The socket attaches firmly at rail ends or notches on American-made operating room tables, providing secure connections with flexible positioning and quick adjustments during procedures.

Constructed from durable stainless steel, the MCM 160 ensures long-lasting performance, corrosion resistance, and minimal maintenance. Its robust design fits all American-made operating room tables, helping surgical teams maintain efficiency and precision throughout procedures.

The MCM 160 is engineered to meet strict healthcare industry standards and provides dependable performance in demanding surgical settings.

FEATURES

- Accepts Round Posts
- Constructed from Stainless Steel
- Long-Lasting Investment
- Compatible with all American-made O.R. Tables
- Attaches at Rail Ends or Notches

SPECIFICATIONS

- Accepts Round Posts up to 11/16" Diameter
- Construction: #304 Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)