

# MCM Stainless Steel Picket Fence Leg Holder | MCM 611



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The MCM 611 Picket Fence Leg Holder is a surgical tool designed to hold a patient's leg in place during operating room procedures. Made from stainless steel, it features a unique picket fence design that secures both legs simultaneously, improving surgical prep speed and patient alignment while maintaining proper leg positioning.

The 5/8" support bar adjusts up to 30" to accommodate various patients and procedures. The leg holder attaches securely to standard OR table rails using a Clark socket (sold separately, MCM 160).

The smooth stainless steel surface resists corrosion and allows for quick sterilization between procedures, supporting strict hygiene protocols. The unit is backed by a 2-year limited warranty.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)