

# MCM Single Loop Leg Prepper with Pad | MCM 157



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The MCM 157 Single Loop Leg Prep Holder provides secure and stable leg positioning during surgical prep. Constructed from durable stainless steel, it attaches easily to standard OR table rails using a Clark socket (sold separately). The height adjusts up to 30 inches, offering flexibility for various procedures and patient sizes.

The leg loop features a cleanable pad that supports both comfort and hygiene in the surgical environment. With a weight capacity of 50 pounds, this holder provides strong, dependable support for safe and efficient leg prepping.

Designed for hospitals, surgical centers, and outpatient facilities, the MCM 157 is American-made and built to meet high standards in surgical equipment design and durability.

## **SPECIFICATIONS**

- Construction Material: Stainless Steel
- Mounting Method: Attaches to standard surgical table side rails via Clark Socket
- Height Adjustment Range: Up to 30 inches
- Support Loop: Includes one cleanable loop pad
- Replacement Pads: Model MCM156 – sold in boxes of 12
- Weight Capacity: 50 lbs (22.7 kg)
- Made in the USA

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)