

MCM Total Knee Replacement (TKR) Support with Pad | MCM 612



Stock Photo - Options May Vary

FACTORY NEW

The MCM 612 TKR Foot Rest provides precise leg support during total knee replacement surgeries. Built from stainless steel and black polyoxymethylene, it delivers strength and stability in the operating room. The radiolucent top allows for clear imaging without repositioning the patient, which is critical for surgical accuracy.

The MCM 612 footrest comes with one black pad and features a 5/8" stainless steel post that fits firmly into any standard Clark socket. It is easy to set up and stays steady, working well with most operating room tables.

Engineered for high-demand surgical environments, the MCM 612 system helps surgical teams improve patient alignment during total knee arthroplasty. Its corrosion-resistant materials ensure long-term use and simple cleaning. Replacement pads are available separately for convenient maintenance between procedures.

SPECIFICATIONS

- Material: Stainless Steel
- Support Bar Diameter: 5/8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)