

MCM Replacement Pad for Total Knee Replacement Support | MCM 613



Stock Photo - Options May Vary

FACTORY NEW

The MCM 613 Total Knee Support Replacement Pad is a precision component designed for use with the MCM612 knee support system during total knee replacement procedures. The pad features a radiolucent top portion that allows for clear imaging without requiring patient repositioning, ensuring optimal visualization during surgery.

This replacement pad is manufactured from black polyoxymethylene material and seamlessly integrates with the MCM612 system to restore full functionality. It is sold as a case of 12 pads, providing a reliable supply for surgical facilities.

The MCM 613 is suitable for orthopedic surgical centers and hospitals requiring consistent, durable components for knee support systems. As an original factory part, it delivers reliable performance for total knee replacement procedures.

SPECIFICATIONS

- Compatible with the MCM612 knee support system (pad only; support system sold separately)
- Made of black polyoxymethylene material
- Sold as case of 12

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)