

Steris-Birkova Continuous Wrap Restraint: PR-3920



Stock Photo - Options May Vary

FACTORY NEW

The Steris-Birkova Continuous Wrap Restraint PR-3920 delivers strong performance and dependable security for surgical and treatment settings. It works with a wide range of tables, including surgical, pain management, IA, and treatment tables, allowing medical teams to streamline patient positioning with greater confidence.

The simple wrap design allows fast adjustments during procedures, helping staff maintain focus and control. The continuous wrap restraint offers flexible sizing to meet different patient and table requirements, with the standard model measuring 80 inches in overall length and supporting secure wrapping without complexity.

Each strap cleans easily and withstands regular use, allowing teams to maintain high standards with less effort. For facilities needing extended reach, 92-inch and 120-inch versions are also available, matching the design of the standard model while providing additional length for broader applications.

FEATURES

- 4" Wide Vinyl
- Velcro Closures
- Wraps Completely Around Table
- Straps Secure Back on Themselves

SPECIFICATIONS

- Construction: 4" Wide Vinyl
- Overall Length: 80"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)