

SchureMed Extra Wide Armboard | 800-0018



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Extra Wide Armboard is a surgical armboard designed for operating room personnel. Nearly twice as wide as the standard model, it accommodates a broader range of patient needs while maintaining all essential features of the standard armboard. The board is lightweight, extremely durable, and radiolucent, making it suitable for use during imaging procedures.

The armboard connects securely to standard U.S. operating room table side rails. A finger-trigger mechanism at the end allows for easy positioning and repositioning. Simply pull the trigger to move the armboard to the desired location, then release it to lock it firmly in place.

The device is designed for single-person installation and requires no maintenance. It is not made with natural rubber latex and has a five-year lifespan under normal use conditions.

FEATURES

- Extra-Wide Armboard
- Connects to O.R. table side rails and locks securely
- To swing the armboard, pull the finger trigger at the end of the armboard and move to the desired position
- Release the trigger, and it locks tight
- Pad Sold Separately (see Options)

SPECIFICATIONS

- Dimensions: 11" x 26.5"L
- Length: 26.5" +/- 0.5"
- Width: 11" +/- 0.5"
- Depth: 2.5" +/- 0.5" (with pad)
- Device Weight: 9 +/- 0.5 lbs
- Compatibility: Fits standard U.S. operating room table side rails
- Weight Capacity: Device supports 500 lb proportional patient load 6'4" tall patient
- The life of the device is five years under normal use

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)