

Steris-Birkova Deluxe Restraint Strap: PR-3400



Stock Photo - Options May Vary

FACTORY NEW

The Steris-Birkova PR-3400 is a deluxe restraint strap designed for surgical and procedure tables. It features a sturdy 4-inch width with two-piece vinyl and Velcro construction that provides a secure hold during procedures and withstands repeated use. The strap cleans easily for efficient maintenance.

The PR-3400 fits a wide range of surgical and procedure tables, offering flexible positioning options and efficient workflow management. It mounts directly to siderails without hooks, which speeds up preparation and reduces clutter around the table.

This reusable restraint strap delivers reliable security and comfort for surgical teams. The simple attachment method and smooth operation integrate easily into daily routines, helping healthcare professionals maintain control with minimal effort.

FEATURES

- Set of two straps
- Reusable
- 4" Wide, 2 Piece Velcro® Closures
- Mounts to Any Style Table with Side Rails
- No Hooks are necessary

SPECIFICATIONS

- Width: 4" Wide
- Length: Extends 30" to 70"
- Velcro® Closures

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)