

# Birkova Leg/Foot Extension: LS-9500



Stock Photo - Options May Vary

## FACTORY NEW

The Birkova Leg/Foot Extension Model LS-9500 is a surgical table accessory constructed from stainless steel and phenolic materials. It extends the operating room table by 10 inches and can be repositioned to function as a vertical foot support. The unit includes a 2-inch thick pad and mounting clamps for secure attachment.

Designed for quick and easy setup, the LS-9500 is compatible with most surgical tables. The dual-function design provides flexibility for various surgical positioning needs while maintaining a compact footprint when not in use.

When ordering, specify the make and model of your surgical table to ensure proper compatibility and fit.

## FEATURES

- Leg/Foot Extension, Model Number LS-9500
- Adds 10 inches to the length of the O.R. table
- Doubles as foot support in a vertical position
- It comes with a 2" pad and mounting clamps
- Quick and easy to setup

## SPECIFICATIONS

- Construction: Stainless Steel and Phenolic
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)