


Maquet 1130 Solid Leg Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Maquet 1130 Solid Leg Standard Series Replacement Pad Set is a 3-piece pad set designed for low-risk patients and short-duration surgical procedures. The set features a solid leg section constructed with medium-firm high-density foam and measures 2 inches thick overall.

The pads are covered with electrically conductive 'Lectrolite vinyl, providing fluid and stain resistance while maintaining X-ray permeability. The set includes waterfall edge construction and is latex-free with anti-microbial properties.

This replacement pad set is compatible with the Maquet Surgical Table Model 1130 and comes complete with a loop and hook kit for secure attachment.

FEATURES

- For Maquet Table Model 1130
- Designed for short-duration procedures
- Solid Leg Section
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-Free

SPECIFICATIONS

- Construction: Medium Firm High-density Foam
- Exterior Cover: Electrically Conductive 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)