


Maquet Alphastar 1132 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Maquet Alphastar 1132 Standard Series Replacement Pad Set is designed to support low-risk patients during short-duration surgical procedures. The pad set features 2-inch thick construction with medium-firm high-density foam and is covered with conductive Lectrolite vinyl for safety and durability.

This complete pad set includes four cushions with loop and hook fastening systems for secure attachment to the surgical table. The waterfall edge construction provides comprehensive patient support and positioning flexibility during procedures.

The pads are fluid and stain resistant, X-ray permeable, and treated with anti-microbial protection. The latex-free design accommodates patients with latex sensitivities.

FEATURES

- Designed for short-duration procedures
- Loop & Hook kit comes complete with the Pad
- 4 Cushion Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-density Foam
- Exterior Cover: Electrically Conductive Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)