


Maquet Alphastar 1132 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Maquet/Alphastar 1132 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features a 2-inch overall thickness constructed with medium-firm high-density foam and 1-inch Visco memory foam, providing optimal support and pressure distribution.

The set includes four cushion pads covered with conductive 4-way stretch Fusion III soft vinyl, designed to accommodate various surgical table configurations. The waterfall edge construction ensures seamless transitions between pad sections for enhanced patient positioning and comfort.

The pads are fluid and stain resistant, X-ray permeable, and treated with anti-microbial protection. Each replacement set comes complete with loop and hook attachment hardware for secure installation on compatible surgical tables.

FEATURES

- Designed for longer-duration procedures
- Loop & Hook kit comes complete with the Pad
- 4 Cushion Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)