

# Maquet Getinge Alphamaxx Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Maquet Getinge Alphamaxx Standard Series Replacement Pad Set is designed to support low-risk patients during short-duration surgical procedures. The set features 2-inch thick padding constructed with medium-firm high-density foam and covered with conductive Lectrolite vinyl for optimal patient comfort and safety.

This complete pad set includes four cushions with waterfall edge construction to prevent fluid pooling and migration. The pads are fluid and stain resistant, X-ray permeable, and treated with anti-microbial protection to maintain a clean surgical environment.

The set is latex-free and includes a hoop and loop kit for secure attachment to the surgical table. The electrically conductive cover ensures proper grounding during electrosurgical procedures.

## FEATURES

- For the Maquet/Getinge Alphamaxx
- Designed for short-duration procedures
- Comes complete with the Pad
- 4 Cushion Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial

## SPECIFICATIONS

- Construction: Medium Firm High-density Foam
- Exterior Cover: Electrically Conductive Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)