

Maquet Getinge Alphamaxx Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Maquet Getinge Alphamaxx Premium Series Replacement Pad Set is designed to enhance patient comfort during longer-duration surgical procedures. The pad set features a 2-inch overall thickness constructed with medium-firm high-density foam and 1-inch Visco memory foam, providing optimal support and pressure distribution.

The pads are covered with conductive 4-way stretch Fusion III soft vinyl that resists fluid and stains while maintaining X-ray permeability. The complete set includes four cushions with waterfall edge construction and a Loop & Hook kit for secure attachment to the surgical table.

The pads are latex-free and treated with anti-microbial protection to maintain a hygienic operating environment. This replacement pad set is compatible with Maquet Getinge Alphamaxx surgical tables.

FEATURES

- Designed for longer-duration procedures
- Loop & Hook kit comes complete with the Pad
- 4 Cushion Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)