

Skytron 6500B Surgical Table Replacement Hand Control



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6500B Surgical Table Replacement Hand Control is a direct replacement for the original hand control on Skytron Model 6500B Elite battery-operated surgical tables. This control unit restores full functionality to tables with non-responsive or malfunctioning hand controls, allowing surgeons and OR staff to resume normal table operation without costly equipment replacement.

The hand control features a 24-pin connector that provides 100% compatibility with Skytron 6500B Elite tables. The unit includes a rail hook for convenient mounting and storage when not in use. Durable construction ensures reliable performance in demanding surgical environments.

Each replacement hand control is backed by a 90-day ship back warranty, providing peace of mind and protection against defects. Before ordering, verify that your table's connector is a 24-pin type to ensure proper compatibility.

FEATURES

- Durable Construction
- 24 Pin Connector
- Rail Hook
- 100% Compatibility
- 90 Day Ship Back Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)