

Skytron 5001 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 5001 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of medium-firm high-density foam combined with 1 inch of Visco memory foam, covered with conductive 4-way stretch Fusion III soft vinyl for durability and performance.

This replacement pad set is compatible with Skytron models 5001, 6001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD. The construction provides reliable support while maintaining the conductive properties necessary for safe surgical table operation.

The pad set includes a complete Velcro kit for installation and is designed to withstand the demands of the operating room environment with fluid and stain resistance.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)