

# Skytron 5001 Surgical Table Replacement Hand Control



*Stock Photo - Options May Vary*

## FACTORY NEW

The Skytron 5001 Surgical Table Replacement Hand Control is a compatible replacement control unit designed for Skytron Model 5001 Elite non-battery operating tables. This hand control restores full functionality to surgical tables with non-responsive or damaged remote controls.

The replacement features a 19-pin connector that ensures proper integration with compatible Skytron tables. A rail hook attachment allows secure mounting during surgical procedures.

The unit is backed by a 90-day ship back warranty, providing assurance of product reliability and manufacturer support.

## FEATURES

- Durable Construction
- 19 Pin Connector
- Rail Hook
- 100% Compatibility
- 90 Day Ship Back Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)