

Skytron 6001 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6001 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2-inch thick medium-firm high-density foam construction covered with electrically conductive vinyl, providing a durable and functional surface for surgical applications.

This replacement pad set is compatible with Skytron models 5001, 6001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD. The complete set includes a Velcro kit for secure attachment to the surgical table.

FEATURES

- Designed for Short-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)