

Skytron 6001 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6001 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. It features a 2-inch thick construction combining medium-firm high-density foam with 1-inch Visco memory foam, covered in conductive 4-way stretch Fusion III soft vinyl.

This pad set is compatible with Skytron models 5001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD. The complete kit includes a Velcro attachment system for secure installation on the surgical table.

The pad construction provides fluid and stain resistance while maintaining X-ray permeability for intraoperative imaging. Anti-microbial properties and latex-free materials support infection control protocols in the operating room.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)