

# Skytron 6500HD Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron 6500HD Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. It features 2-inch thick construction with medium-firm high-density foam covered in electrically conductive vinyl. This pad set is compatible with Skytron models 5001, 6001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD.

The pad set comes complete with a Hook & Loop kit for secure attachment. The electrically conductive vinyl exterior provides fluid and stain resistance while maintaining X-ray permeability for imaging compatibility. The waterfall edge construction ensures proper fit and patient positioning on the surgical table.

The pad is latex-free and features anti-microbial properties to support infection control protocols in the operating room.

## FEATURES

- Designed for Short-Duration Procedures
- It comes Complete with a Pad and Velcro kit
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Hook & Loop Kit Included

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)