

Skytron 6500 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6500 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. It features a 2-inch thick construction combining medium-firm high-density foam with 1-inch Visco memory foam, covered in conductive 4-way stretch Fusion III soft vinyl.

This pad set is compatible with Skytron models 5001, 6001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD. The Velcro kit comes complete with the pad for easy installation and secure positioning on the surgical table.

The pad surface is fluid and stain resistant, X-ray permeable, and treated with anti-microbial protection. It is latex-free and features waterfall edge construction for enhanced durability and patient safety.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)