

# Skytron 6600 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron 6600 Standard Series Replacement Pad Set is designed for low-risk patients undergoing short-duration procedures. The pad features 2 inches of medium-firm high-density foam construction covered with electrically conductive vinyl, providing reliable support and protection during surgical use.

The pad set includes waterfall edge construction and is fluid and stain resistant, making it easy to maintain and clean between procedures. The electrically conductive cover ensures compatibility with electrosurgical equipment while the X-ray permeable design allows for imaging without pad removal.

The set comes complete with a Hook and Loop attachment kit for secure positioning on the surgical table. The pad is latex-free and treated with anti-microbial protection to support infection control protocols.

## FEATURES

- Designed for Short-Duration Procedures
- Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Includes Hook & Loop Kit.

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)