

Skytron 6600 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6600 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features a 2-inch overall thickness combining medium-firm high-density foam with 1-inch Visco memory foam for optimal support and pressure distribution. The conductive 4-way stretch Fusion III soft vinyl cover provides durability and ease of cleaning while maintaining compatibility with surgical equipment.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)