

Skytron 6700B Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6700B Standard Series Replacement Pad Set is designed for low-risk patients and short-duration surgical procedures. This pad set is compatible with Skytron models 6700B and 6701 surgical tables. The set features 2-inch thick construction with medium-firm high-density foam and an electrically conductive vinyl cover.

The pad is constructed to provide essential support during routine procedures while maintaining durability and ease of cleaning. The waterfall edge construction and fluid-resistant properties help protect the underlying table structure. The set includes a complete hook and loop fastening kit for secure installation.

For patients requiring extended procedure times or additional comfort needs, alternative premium pad options are available. This standard pad set is latex-free and includes anti-microbial properties to support infection control protocols in the operating room.

FEATURES

- Designed for Short-Duration Procedures
- Velcro Kit comes Complete with Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex-free
- Hook & Loop Kit included

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)