

# Skytron 6700B Premium Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron 6700B Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. This pad set is compatible with Skytron models 6700B and 6701 surgical tables. The 2-inch thick pad combines medium-firm high-density foam with 1-inch visco memory foam for optimal support and pressure distribution.

The pad surface is covered with conductive 4-way stretch Fusion III soft vinyl, providing durability and ease of cleaning. The waterfall edge construction ensures secure positioning and prevents fluid migration. The cover material is fluid and stain resistant, making it suitable for demanding operating room environments.

The pad set includes a complete Velcro kit for installation and features anti-microbial properties and X-ray permeability. The construction is latex-free, accommodating patients with latex sensitivities.

## FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)