

# Skytron 6701 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron 6701 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration surgical procedures. This pad set is compatible with Skytron models 6700B and 6701 surgical tables. The set features 2-inch thick construction with medium-firm high-density foam and an electrically conductive vinyl cover.

The pad provides essential protection and comfort during standard surgical cases. Its waterfall edge construction and fluid-resistant properties help maintain a clean surgical environment. The latex-free design accommodates patients with sensitivities.

The complete pad set includes hook and loop fasteners for secure attachment to the surgical table. The X-ray permeable construction allows for imaging procedures without pad removal.

## FEATURES

- Latex-Free
- Designed for Short-Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Hook & Loop included
- Velcro Kit Complete with Pad

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)