

Skytron 6701 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6701 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features a 2-inch overall thickness constructed with medium-firm high-density foam and 1-inch Visco memory foam, covered with conductive 4-way stretch Fusion III soft vinyl. This replacement pad set is compatible with Skytron models 6700B and 6701 surgical tables.

The premium construction provides enhanced support and comfort for longer-duration cases while maintaining durability and ease of cleaning. The conductive cover and fluid-resistant properties make it suitable for demanding operating room environments.

The complete pad set includes a Velcro kit for secure attachment to the surgical table, ensuring stability during procedures.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)