

# Skytron Elite 6002 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron Elite 6002 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2 inches of medium-firm high-density foam construction covered with electrically conductive vinyl for reliable performance in the operating room.

This replacement pad set includes a complete waterfall edge construction and comes with a hook and loop kit for secure installation. The design prioritizes essential functionality and durability for standard surgical applications.

## FEATURES

- Designed for Short-Duration Procedures
- Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free
- Hook & Loop Kit is included

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)