

Skytron 6002 Surgical Table Replacement Hand Control



Stock Photo - Options May Vary

FACTORY NEW

Replacement hand control for Skytron Model 6002 Elite surgical tables. This control unit restores full functionality to non-battery-operated tables when the original hand control is damaged or malfunctioning. The unit features a 19-pin connector that directly interfaces with compatible Skytron 6002 Elite systems.

Designed for durability and reliable performance in the operating room environment. The hand control includes a rail hook for convenient mounting and storage during procedures. Installation is straightforward for tables equipped with the matching 19-pin connector interface.

Each replacement unit comes with a 90-day ship back warranty covering defects and operational issues. Verify that your table's connector is a 19-pin type before ordering to ensure compatibility.

FEATURES

- Durable Construction
- 19 Pin Connector
- Rail Hook
- 100% Compatibility
- 90 Day Ship Back Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)