

Skytron Elite 6002 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 6002 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of total thickness, combining medium-firm high-density foam with 1 inch of visco memory foam for optimal support and pressure distribution.

The pad is covered with conductive 4-way stretch Fusion III soft vinyl that provides durability and ease of cleaning. The waterfall edge construction and fluid-resistant properties make it suitable for demanding operating room environments while maintaining infection control standards.

This replacement pad set comes complete with a velcro kit for secure attachment to the surgical table. The design is latex-free and incorporates anti-microbial properties to support a hygienic surgical environment.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)