

Skytron Elite 3600B Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 3600B Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2 inches of medium-firm high-density foam construction covered with electrically conductive vinyl, providing a balance of support and comfort for routine surgical applications.

The pad set includes waterfall edge construction and is fluid and stain resistant, making it easy to clean and maintain. It is X-ray permeable, anti-microbial, and latex free, meeting standard infection control and safety requirements for operating room use.

FEATURES

- Designed for Short-Duration Procedures
- Replacement Pad Comes Complete with Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)